





Vysis ALK Break Apart FISH Probe Kit
VP 2000 Reference Guide



CHOOSE TRANSFORMATION

VP 2000 PRETREATMENT PROTOCOL AND REAGENT MAP

Step	Description	Basin	Mode	Reagent	Time (m:s)	Temp (°C)
1	Hemo-De	4	—	Hemo-De	5 minutes	Ambient
2	Hemo-De	5	—	Hemo-De	5 minutes	Ambient
3	Hemo-De	6	—	Hemo-De	5 minutes	Ambient
4	100% Ethanol	7	—	100% Ethanol	1 minute	Ambient
5	100% Ethanol	8	—	100% Ethanol	1 minute	Ambient
6	Pretreatment SSC Solution	1	—	Pretreatment SSC Solution	25 ± 10 minutes	80°C
7	Purified Water	9	—	Purified Water	3 minutes	Ambient
8	Protease Solution	3	—	Protease Solution	15 ± 5 minutes	37°C
9	Purified Water	10	—	Purified Water	3 minutes	Ambient
10	70% Ethanol	12	—	70% Ethanol	1 minute	Ambient
11	85% Ethanol	13	—	85% Ethanol	1 minute	Ambient
12	100% Ethanol	14	—	100% Ethanol	1 minute	Ambient
13	Air Dry	Drying Station	—	—	—	Ambient



VP 2000 SLIDE WASHING PROTOCOL AND REAGENT MAP

Step	Description	Basin	Mode	Reagent	Time (m:s)	Temp (°C)
1	Wash Buffer II	2	—	Wash Buffer II	2 minutes	74°C
2	Air Dry	Drying Station	—	—	30 seconds	Ambient

